



Intelligent Flow Peristaltic Pump LIFP-A1 series (One Channel)

Intelligent Flow Peristaltic Pump LIFP-A1 series

Intelligent Flow Peristaltic Pump LIFP-A1 series are equipped with a touch screen controller that allows monitoring of parameters and change of speed. Equipped with advanced features such as flow rate calibration and anti-drip function. Easy dispense mode to dispense the recorded volume by one touch DISPENSE key. Intelligent cooling fan control helps in minimizing the work noise

Features

- LCD display for monitoring parameters
- Anti-drip function ensuring dispense accuracy
- Speed resolution : 0.1 rpm
- Adjustable rotating speed
- Intelligent temperature control minimizing working noise
- Reversible operation allowing pumping in both directions
- Multi channels with interchangeable pump heads for varied applications
- Easy liquid auto dispensation by time or volume
- Internal double layer isolation
- Dust and moisture proof circuit board
- Optional remote infrared control
- Robust and protected with ABS plastic housing
- Low noise operation

Intelligent Flow Peristaltic Pump LIFP-A1 series

Specification

Model	LIFP-A10	LIFP-A11
Flow Range	0.006 ~ 1600 mL/min	0.006 ~ 2900 mL/min
Flow Accuracy	0.5 %	0.2 %
Speed Range	0.1 ~ 350 rpm	0.1 ~ 600 rpm
No of channel	1	1
Speed Resolution	0.1 rpm	0.1 rpm
Ambient Temperature	0 ~ 40°C	0 ~ 40°C
Relative Humidity	< 80 %	< 80 %
IP grade	IP 31	IP 31
Power supply	AC 220 V \pm 10 %, 50 Hz / 60 Hz	AC 220 V \pm 10 %, 50 Hz / 60 Hz
Power Consumption	< 40 W	< 50 W
Dimensions (L x W x H)	257 x 198 x 182 mm	257 x 198 x 182 mm
Weight	5.3 kg	5.5 kg

Pump Heads and Tubing

Model No	Suitable Pump Head	Channel	Tubing Size (mm)	Flow rate (mL/min) per channel
LIFP-A10	AB1	1	Wall : 0.8 ~ 1	0.006 - 990
	AB2	1	Wall : 0.8 ~ 1	0.16 - 990
	AB3	1	Wall : 0.8 ~ 1	0.006 - 1300
	AB4	1	Wall : 0.8 ~ 1	0.16 - 1600
LIFP-A11	BA1	1	Wall : 0.8 ~ 1	0.006 - 1700
	BA2	1	Wall : 0.8 ~ 1	0.16 - 1700
	BA3	1	Wall : 0.8 ~ 1	0.006 - 2300
	BA4	1	Wall : 0.8 ~ 1	0.16 - 2900