



-25°C Chest
Freezers

www.labtron.com

info@Labtron.com

-25°C Chest Freezers

-25°C Chest Freezers LCF-A1 Series

Chest type low temperature freezer comes in top door key locking facility. The -25°C low temperature freezer has CFC free compressor and quick freezing switch for fast cooling. They are equipped with German EBM fan motor and R134a as Refrigerant. These freezers comes in 110L, 196L, 226L, 358L, 508L capacities and are ISO9001, ISO14001, ISO13485 certified.

Features


- Contains one basket made of steel wire for storage and four casters for easy handling
- Outer body made up of stainless steel with aluminum panel inside the freezer
- CFC free compressor with quick freezing switch for fast cooling
- Equipped with German EBM fan motor and R134a as Refrigerant
- It has microprocessor-based temperature controller with range -10°C to -25°C
- User friendly with digital temperature display
- Perfect audio/visual alarm for high-low temperature and system failure

Applications

The -25°C chest type low temperature freezers are designed to store large size items, frozen eatables, reagents, chemicals & injectable along with other goods that require freezing storage temperature.

Specifications

Model	LCF-A10	LCF-A11	LCF-A12	LCF-A13	LCF-A14
Capacity	110L	196L	226L	358L	508L
External Dimension (W×D×H)	549×549×845mm	556×1056×937mm	556×1206×937mm	730×1204×968mm	730×1554×968mm
Internal Dimension (W×D×H)	410×410×654mm	430×930×480mm	430×1080×480mm	530×1080×625mm	530×1400×685mm
Input Power	145W	180W	207W	320W	375W
Weight (Nt/Gt in Kg)	30/40	50/60	55/65	80/90	100/110
Power Supply	220V/50Hz				



Labtron Equipment Ltd.
Sentinel House, Ancells Business Park
Harvest Crescent, Fleet, GU51 2UZ, UK
Tel: 01252 413773
www.labtron.com • info@Labtron.com