



LABTRON



-25°C Upright
Freezers

www.labtron.com

info@Labtron.com

-25°C Upright Freezer

-25°C Upright Freezer LUF-A2 Series

The -25°C upright type low temperature freezer comes in safety lock on door with handle. It has CFC free compressor and quick freezing switch for fast cooling. These upright freezers comes with 4 caster unit for easy handling and 7-12 units drawers for more storage. The -25°C upright freezer with 270 Liters and 450 Liters capacity are designed to efficiently store eatables, reagents, chemicals & injectables along with other goods that require freezing storage temperature. They are ISO9001, ISO14001, ISO13485 certified freezers.

Features

- Equipped with microprocessor temperature controller with temperature range -10°C to -25°C
- High or low temperature alarm and system failure alarm with audible/visual alarm function
- Made up of ABS material from inside with stainless steel outer body
- Comes with 4 caster units for easy handling and 7-12 units drawers for more storage
- Fully automatic, LCD display
- CFC free compressor with quick freezing switch for fast cooling
- Equipped with German EBM fan motor and R507 Refrigerant

Applications

The -25°C upright type low temperature freezers are suitable for use in food industries, blood banks, medical laboratories and research institutes.

Specifications

Model	LUF-A20	LUF-A21
Capacity	270L	450L
External Dimension (W×D×H)	700×600×1750mm	655×780×1845mm
Internal Dimension (W×D×H)	440×400×1490mm	530×626×1356mm
Input Power	370W	450W
Weight (Nt/Gt in Kg)	65/70	120/140
Power Supply	220V/50Hz	



Labtron Equipment Ltd.

**Sentinel House, Ancells Business Park
Harvest Crescent, Fleet, GU51 2UZ, UK**

Tel: 01252 413773

www.labtron.com • info@Labtron.com