

Chemical

STORAGE CABINET



Chemical Storage Cabinet LCSC-A1 Series

The chemical storage cabinets are controlled by microprocessor and allows automatic discharge of toxic fumes and pungent gases via the duct and filter as clean air. These chemical storage cabinets comes with dark compartment for storing light sensitive chemicals. Powder coated steel exterior makes these chemical storage cabinets corrosion resistant. They are equipped with casters for easy installation. They are safe, reliable and consists of recovery system in case of sudden power cut.

Features

- Different set of chambers & adjustable shelves
- LCD digital display and LED lamp
- Recovery system in case of sudden power cut
- Separate door and lock system for each compartment
- Airtight door prevents leakage of toxic gases
- Noiseless blower motor
- Auto temperature control system



Application

Chemical Storage Cabinets are used for storing corrosive and toxic chemicals in chemical laboratories, pharmaceutical industry, biochemistry, biotechnology, environmental science.



Specifications

Model	LCSC-A10	LCSC-A11	LCSC-A12	LCSC-A13	LCSC-A14
No. of Chambers	4 + 2 dark	4 + 1 dark + 1 low temperature	2 + 1 low temperature	2 + 1 dark	2
Controller	Digital LCD Display	In-built CPU/Auto temperature control system		In-built CPU/Auto temperature control system	
Overall Dimension	1240x760x2200 mm	1240x760x2200 mm	820x760x2250 mm	820x760x2250 mm	900x760x2050 mm
Power	265W	390W	140W	80W	75W
No. of Shelves	14		7		6
Blower Motor	Sirocco Fan for noiseless operation				
Filter	Carbon Filter				
Lamp	LED (High luminance)				
Interior material	ABS injection molding				
Exterior material	Steel Plate with Powder Coating				
Timer	10min on / 10min off				
Power Supply	AC 230V, 50/60Hz				



Labtron Equipment Ltd

Basing View, Basingstoke, Hampshire RG21 4RG, UK
01256 610 110 | info@Labtron.com
www.Labtron.com