



# Cube Ice Makers (A Series)



## Cube Ice Makers-A series

Cube Ice maker series comes with insulated ice storage bin and door which minimizes air condensation and water formation on equipment surface. It produces steady and similar shape ice cubes with high transparency. It serves safe, hygienic, odorless and best tasting cube ices which can be used for all types of functions. The machine is configured using high quality material and advanced technology at par with ISO and CE standards. It has integrated design for maximum energy saving and water conservation.

---

### Features

- ✔ Constructed with stainless steel and rigid plastic cabinets and other high-quality components
  - ✔ Auto indicator for indicating level of ice and water
  - ✔ Ice shape is steady and similar with high transparency
  - ✔ R404A as refrigerant, high efficiency CFC free compressor
  - ✔ Integrated design for maximum energy saving and water conservation
  - ✔ User-friendly electronic controls, overload protection
  - ✔ Produces edible ice using potable water
  - ✔ Testing conditions: Ambient temp. 20°C, Inlet water temp. 10°C
- 

### Applications

Cube ice makers are widely used in restaurants and other food service operations for cooling water and beverages. They are ideal for bars to enhance the appearance of alcoholic drinks due to clarity of ice. Often used in laboratories, hotels, marine fisheries, chemical industries, pharmaceutical industries etc.

# CUBE ICE MAKERS

## Specification

Model	LCIM-A10	LCIM-A11	LCIM-A20	LCIM-A21	LCIM-A22
Capacity	40kgs/24hrs	50kgs/24hrs	60kgs/24hrs	75kgs/24hrs	100kgs/24hrs
Storage Bin Capacity	20kgs	10kgs	38kgs	38kgs	10kgs
Dimension (WxDxH)mm	635x570x830mm		765x740x900mm		
Cooling Mode	Air cooling				
Cooling Fan	EBM(Germany)				
Dry Filter	Danfoss(Denmark)				
Power Consumption	420W	440W	460W	480W	
Standard Voltage	1P/220V/50Hz				
Weight	44kgs		61kgs		
Exterior Finish	304/2B Stainless steel				
Compressor Brand	Aspera (Slovenia)				
Refrigeration System	R404A				
Certification	CE				

Model	LCIM-A23	LCIM-A24	LCIM-A30	LCIM-A31	LCIM-A32
Capacity	125kgs/24hrs	145kgs/24hrs	160kgs/24hrs	195kgs/24hrs	230kgs/24hrs
Storage Bin Capacity	58kgs	58kgs	82Kgs	120kgs	236kgs
Dimension (WxDxH)mm	765x740x900mm	635x570x830mm		765x740x900mm	
Cooling Mode	Air cooling				
Cooling Fan	EBM(Germany)				
Dry Filter	Danfoss (Denmark)				

# CUBE ICE MAKERS

## Specification

Power Consumption	730W	750W	860W	950W	1052W
Standard Voltage	1P/220V/50Hz				
Weight	61kgs		83kgs	87kgs	-
Exterior Finish	304/2B Stainless steel				
Compressor Brand	Aspera (Slovenia)				
Refrigeration System	R404A				
Certification	CE				

Model	LCIM-A33	LCIM-A34	LCIM-A35	LCIM-A36	LCIM-A37
Capacity	295kgs/24hrs	340kgs/24hrs	454kgs/24hrs	545kgs/24hrs	636kgs/24hrs
Storage Bin Capacity	236kgs		280kgs		
Dimension (WxDxH)mm	765x740x900mm		760x800 x1750mm (1.1CBM)	844x705x620mm; 840x790x1055mm (1.1 CBM)	
Cooling Mode	Air cooling				
Cooling Fan	EBM (Germany)				
Dry Filter	Danfoss (Denmark)				
Power Consumption	1150W	1250W	1600W	1980W	2600W
Standard Voltage	1P/220V/50Hz				
Weight	105kgs	107kgs	120kgs	122kgs	128kgs
Exterior Finish	304/2B Stainless steel				
Compressor Brand	Aspera (Slovenia)				
Refrigeration System	R404A				
Certification	CE				



**Labtron Equipment Ltd**

**Sentinel House, Ancells Business Park,**

**Harvest Crescent, Fleet**

**GU51 2UZ, UK**

**Telephone: 01252 413773**

**[www.labtron.com](http://www.labtron.com) • [info@Labtron.com](mailto:info@Labtron.com)**