



PORTABLE AUTOCLAVE

LPOA-A10

www.labtron.com

info@labtron.com

Portable Autoclave LPOA-A10

Portable autoclave LPOA-A10 is a portable N class unit with capacity upto 12 L and a sterilizing temperature of 126 °C. It has an effectively sterilization which can inactivate even HIV, HBV, mad cow virus and bacillus. Equipped with safety features like automatic shutoff preventing overheating of chamber and auto protection against raised temperature and pressure.

Features

- Working temperature : 126 °C
- Equipped with interior circulate system
- Self-inflating leak proof chamber
- Safe door lock system ensuring safety / protection during sterilization
- Automatic power cut off with alarm with low water levels indication
- Double scale indicating pressure gauge
- Electric / LPG heated

Application

Used in industries to sterilize equipment, bacteriological and research laboratories, medical, clinical and pharmaceutical industries

Specifications

Model No	LPOA-A10
Volume of Chamber	12 L
Working temperature	126 °C
Heating mode	Electric / LPG heating
Working pressure	0.14 - 0.16 MPa
Max safety pressure	0.165 MPa
Heat average	$\leq \pm 1^{\circ}\text{C}$
Temp adjustment	105 - 126 °C
Power	AC 220 V / 50 Hz
Overall Dimension	360 × 360 × 420 mm
Net Weight	11 Kg
Gross Weight	12 Kg