


Rotary Vane VACUUM PUMP



ROTARY VANE VACUUM PUMP LRVP-A SERIES

Our model of Rotary Vane Vacuum Pump LRVP-A series is designed with unique lubrication system, built-in oil pump and have a large oil volume. In addition to providing smooth, pulse-free air flow, they are compact and simple to install and operate, have low starting and running torque requirements, produce minimal noise or vibrations and are capable of long, service-free life.

Features

- ▶ Compact, reliable and powerful
- ▶ Space-saving and ideal for systems integration
- ▶ High pumping speeds with small dimensions
- ▶ Gas ballast valve to pump little water vapour
- ▶ Oil-anti-suck back device for convenient starting
- ▶ Efficient air cooling system
- ▶ Eco-friendly exhaust gas free of oil mists
- ▶ Quiet, low-vibration operation
- ▶ Easy to maintain

Applications

This product is used for vacuum packaging, drying, distillation, plastic moulding etc..



Specification

Model	LRVP-A10	LRVP-A11	LRVP-A12	LRVP-A13	LRVP-A14
Pumping Speed	166.67L/min	333.33L/min	666.67 L/min	1050L/min	1666.67L/min
Ultimate Pressure	80Pa	80Pa	50 Pa	50Pa	50Pa
Inlet and Outlet	G1/2inch	G1/2inch	G11/4inch	G11/4inch	G11/4inch
Weight	16kg	18kg	40/58kg	45kg	115kg
Rotating Speed	1400rpm	2800rpm	1400 rpm	1400rpm	1420rpm
Oil Capacity	0.5L	0.5L	2 L	2L	4L
Dimension	410x280x210mm		630x290x260 mm	660x300x270 mm	770x480x300 mm
Power	370/550W	750/900W	1100W	1500W	2200W
Power Supply	220V, 50Hz	220V, 50Hz	220 V, 50 Hz	220V, 50Hz	220V, 50Hz



LABTRON EQUIPMENT LTD
Basing View, Basingstoke, Hampshire RG21 4RG
UK

 01256 610 110 |  info@Labtron.com |  www.Labtron.com