



TABLE TOP AUTOCLAVE

LTTA-A1 SERIES

Table top Autoclave LTTA-A1 Series

Table top autoclave LTTA-A1 series is robust and compact, space efficient unit with 20 L of capacity and a sterilizing temperature of 134 °C. It is equipped with 4 sterilizing baskets along with safety features like automatic shutoff preventing overheating of chamber, automatic cold air discharge and auto protection against raised temperature and pressure making a reliable and user friendly device.

Features

- ▶ Digital LCD display to monitor working status and parameters
- ▶ Working temperature : 134 °C
- ▶ Automated drying function
- ▶ Door lock feature : Door cannot be opened till the pressure inside the chamber doesn't fall to 0.027 MPa
- ▶ Safe door lock system ensuring safety / protection during sterilization
- ▶ Automatic power cut off with alarm with low water levels indication
- ▶ Stainless steel sterilizing baskets
- ▶ Automatic cold air discharge from system chamber ensures dry and clean environment
- ▶ Auto exhaustion of cool air

Application

Used to sterilize culture media, instruments, dressings, intravenous equipments, rubber closure of glass equipments used in aseptic techniques and numerous objects that can withstand high temperature and pressure

Specification

Model	LTA-A10	LTA-A11	LTA-A12	LTA-A13
Capacity	20L	24L	35L	50L
Sterilization Temperature	134°C			
Sterilization Pressure	0.22Mpa			
Timer Range	0-60 Minutes			
Heat Average	±1°C			
Chamber Dimension (mm)(Φ×D)	250×420	280×520	300×500	340×550
Net Weight (Kgs)	38	44	80	93
Gross Weight (Kgs)	43	50	95	112
Power Supply	1.5KW/AC220V 50Hz	1.5KW/AC220V 50Hz	2.7KW/AC220V 50Hz	2.7KW/AC220V 50Hz