



# Vacuum Concentrator Centrifuge LVC-A1 series

[www.labtron.com](http://www.labtron.com)  
[info@Labtron.com](mailto:info@Labtron.com)



## Vacuum Concentrator Centrifuge

Vacuum concentrator centrifuge is specifically designed for centrifugation, vacuumizing and heating to concentrate the samples. It has exquisite design, LCD touch screen monitor with micro-computer control system and direct operation interface. Wide range of rotors are available which are compatible with this vacuum centrifuge.

---

### Features

- ▶ Equipped with TFT true colour wide screen touch monitor
- ▶ Aluminum-alloy made diving room
- ▶ Maintenance-free and non-contact drive and rotation system
- ▶ Used in combination with cold-trap
- ▶ Transparent covering for easy operation

---


### Applications

Vacuum Concentrator centrifuge used for removing solvents from samples finds application in molecular biology, genetics and pharmaceuticals.

---

### Specifications

Model No.	LVC-A10	LVC-A11
Rotor capacity	1.5ml×62	10ml×12, 50ml×6
Speed	1350rpm	1500rpm
Temp. Range	Room temp. to 55°C	Room temp. to 60°C
Condensation Time	0-99h59min	
Optional	Vacuum Pump	
Dimension(L×W×H)	500×490×280mm	



Labtron Equipment Ltd.  
Sentinel House, Ancells Business Park  
Harvest Crescent, Fleet, GU51 2UZ, UK  
Tel: 01252 413773  
[www.labtron.com](http://www.labtron.com) • [info@Labtron.com](mailto:info@Labtron.com)